

Product Description:

Gus Gear's LOCK 3000 is a wearable tube and catheter securement system. The device uses an adjustable strap system to secure the tube/catheter by absorbing the force of pulls. The LOCK 3000 is designed for feeding tubes and peritoneal dialysis (PD) catheters and can be adjusted to fit most tubing.

Expected Clinical Benefits:

- Reduced risk of tube or catheter dislodgement and accidental removal
- Decreased tube or catheter movement
- Reduction or elimination of need for adhesive securement products
- Reduction in risk of line entanglement

Intended Use:

The LOCK 3000 is intended for comprehensive tube and catheter securement in:

- Feeding tubes (gastrostomy tubes, jejunostomy tubes, gastrostomy/jejunostomy tubes)
- Peritoneal dialysis (PD) catheters
- Cecostomy tubes

Warnings and Precautions:**Warnings:**

- Intended for **Single Patient Use only**. May be re-used in the same patient for an extended time.
- This product contains nylon, spandex, polyester, polyethylene, polyurethane, medical-grade silicone and calcium carbonate. Do not use in patients with known or suspected allergy to any of these components.
- Do not use if Dual Safe Strip securement strap is missing, damaged or not properly securing the tube or catheter in place for any reason.
- Contact your doctor immediately if the tube or catheter is damaged or dislodged.
- Do not use if patient's body size is too small or too large to fit appropriately into the standard size range of the LOCK 3000.

Precautions:

- Read and follow all instructions prior to use, including care instructions.
- LOCK 3000 must be properly fitted and used according to the instructions to secure tubing/catheter.
- DO NOT remove Dual Safe Strip while tube or catheter is attached. To remove Dual Safe Strip, first unfasten Dual Safe Strip strap(s), releasing tubing/catheter.
- DO NOT place Dual Safe Strip edge in contact with tube/catheter where tube/catheter enters the LOCK 3000 or damage may occur.
- DO NOT use hot water to wash. DO NOT place in a dryer. (Please refer to care instructions.)
- LOCK 3000 should be replaced if hook and loop closures (Velcro) no longer adequately secure adjustment straps.
- Device is NOT sterile.

Possible Complications:

The potential exists for complications including the following:

- If the product is not used within guidelines provided, damage to the catheter or tubing can occur.
- Contact your doctor immediately if the tube or catheter is damaged or dislodged.

MRI Information:

- This product does not contain any metal and is composed entirely of materials that are electrically non-conductive, non-metallic, and non-magnetic.
- This product is considered **MR Safe** but has not been tested for safety with MRI.

Preparation for Use:






- Examine the LOCK 3000 carefully to determine that all parts are present and in working order.
- Do not use if product is damaged or missing parts.
- Product can be re-used for the same patient after washing. (Refer to care instructions.)

Application & Removal:

- Refer to detailed instructions.

Product Tracing:

- Product can be traced based upon a Lot number (Batch Code) located on the Dual Safe Strip adjustable strap.

	Single Patient, Multi-Use
	Refer to IFU for more information
	Not made with Natural Rubber Latex
	Batch Code
	MR Safe

Gus Gear and Break Free! Are trademarks of Gus Gear Inc.
 Copyright 2023 Gus Gear Inc. All rights reserved.

Contact:

Gus Gear Inc
 PO Box 607
 Valencia, PA 16059 USA
 Customer Support Line: 866-822-4GUS (4487)
 Web: <https://gusgear.com/>